



PRIDE®
MOBILITY
LIVE YOUR BEST®

Tel: +31 (0)85-0408630 . Fax: +31 (0)85-0430039
pridemobility.eu . info@prideeurope.nl

Ervaar de schoonheid van lichtgewicht reizen met de i-Go®+. Deze compacte elektrische rolstoel kan in slechts een paar eenvoudige stappen worden opgevouwen voor moeiteloos transport. De i-Go®+ maakt reizen een fluitje van een cent dankzij de handige opbergruimte.



Comfortabele hoofdsteen



Swing away controller



Opbergruimte onder het zadel

Specificaties

Maximale belasting	113 kg
Snelheid	Tot 5,8 km/u
Maximale actieradius	Tot 15 km
Draaicirkel	616 mm
Bodemvrijheid	79 mm
Lengte	1016 mm
Breedte	597 mm
Aandrijf wielen	12" luchtbanden
Zwenkwieken	8" massief
Aandrijving	Twee motoren, achterwielaandrijving
Gewicht zonder accu's	24 kg
Accu's	Type: Lithium ion Capaciteit: 18.2 Ah
Meer?	Download het specificaties blad op de website

Kenmerken

- Een compact en lichtgewicht ontwerp maakt reizen een fluitje van een cent
- Zak op de achterbank voor extra opbergruimte onderweg
- In slechts een paar eenvoudige stappen op te vouwen voor moeiteloos transport
- Groot voetplatform voor het comfort van de gebruiker
- Zitting van speciaal schuim voor extra comfort

1. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
3. Battery weight may vary based on battery manufacturer.
4. Includes standard seat weight.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

Beschikbare kleuren



Zilver